

## CLIA Uromodulin

**EAN-Code:** 8595635311482  
**Katalognummer:** CL-UMOD050  
**Verpackungsgröße:** 50 tests  
**Lagerung:** 2-8 °C  
**Hersteller:** TestLine Clinical Diagnostics s.r.o.



The CLIA immunoassay is based on the principle of magnetic particle separation, and the detection system utilizes flash chemiluminescence with the help of acridinium ester, and the intensity of the light flash is measured in relative light units (RLU).

### General information about the kit:

- Each individual test contains all the necessary reagents (ready-to-use) for conducting the examination.
- Calibration is included in the package.
- Results are evaluated using ng/ml.

For more information, please refer to BiovendorGroup [CLIA webpage](#).